

# ONEcleaner Enzyme Detergent Cleaner Lubricant

Enzyme Surgical Instrument Detergent Lubricant Cleaners **clean faster**.  
The **ONEcleaner** Enzymatic Surgical Instrument Detergents **cut costs**.  
Lubricating Enzyme Surgical Instrument Detergents clean **residue free**.

The neutral pH **ONEcleaner** delivers the four enzymes necessary to break down all forms of surgical proteinaceous bioburden with the surfactant detergents to clean the surface.

- ) Enzyme Cleaners (protease, amylase, carbohydrase, and lipase)
- ) Detergent Cleaners (surfactant cleaners)
- ) Lubricants (water soluble and free rinsing)
- ) Surface Conditioners (strengthen the chromium passive layer of stainless steel)

The **cost cutting** Enzyme Detergent Surgical Instrument Lubricating Cleaners benefits include:

- ) Combined enzyme, detergent, lubricant cleaning treatments improves turnaround times
- ) Replacing multiple products reduces inventory and purchasing costs
- ) High concentration, significantly reduces the amount of product being used
- ) Residue cleaning renders improved instrument performance and longer usable lifespan

The **ONEcleaner** is highly concentrated and delivers an average of 400% more cleaning power. One gallon can deliver 512 gallons (1938 liters) of cleaning power. The dilution rate is .25 to .5 ounces per gallon diluent (1.9 ml liters to 3.8 ml per liter) depending on the level of bioburden, encrustation and stains.

0225-32 easy FOAM-it Pre Wash Cleaner, 12 ea. 32 oz. per case  
0225-1 ONEcleaner, 4 gallons with 1 pump per case  
0225-5 ONEcleaner, 5 gallons  
0225-15 ONEcleaner, 15 gallons



Customer Service: 1-800- 341-7048 email: [CustomerService@mccmed.com](mailto:CustomerService@mccmed.com)  
Technical Support: 1-509-747-5027, email: [Temple@ergo-Logistics.com](mailto:Temple@ergo-Logistics.com)  
Website: [www.all-in-one-surgical-instrument-cleaner.com](http://www.all-in-one-surgical-instrument-cleaner.com)